

MEDIVATORS™
ENDO GATOR™
Irrigation Pump



Automate Irrigation and Reduce Procedure Time

MEDIVATORS™ ENDOGATOR™ Irrigation System eliminates manual syringing. Irrigation is provided via footswitch, reducing procedural time and increasing visibility.

Pump Specifications

Electrical

Voltage/Current Requirement:
120V/240V~50-60Hz

Pump

Flow Rate: 0-550 ml/min
Pressure: 0-48 psi

Ideal Flow Settings

Auxiliary Water Channel: 40-50
Biopsy Channel: 70 to Max

Physical Dimensions

Height: 152 mm Depth: 203 mm
Width: 178 mm Weight: 2.3 kg

MEDIVATORS™ ENDOGATOR™ Accessories

Safety












- Reduced risk of cross contamination
- Endoscope connectors also available in single patient version

Compatibility

- Works with all sterile water bottles
- Connects to all major brands of endoscopes

Efficiency

- Minimizes reprocessing requirements
- Creates consistent practices
- Cost effective and completely disposable
- Endoscope connectors for all methods of irrigation (Biopsy, Jet)

	Tubing	Biopsy Channel Irrigation	Jet Channel Irrigation	
OLYMPUS® GI Endoscopes	200230U 	100303 	100115 	100241 
PENTAX® GI Endoscopes		100136 	100116 	100242 
STORZ FUJINON® GI Endoscopes	200230FU 	100303 	Backflow Valve 100126 	Direct Luer Connection on Tubing 

ENDO GATOR™ and MEDIVATORS™ are trademarks of Medivators Inc.
PENTAX® is a registered trademark of Hoya Corporation.
OLYMPUS® is a registered trademark of OLYMPUS Corporation.
FUJINON® is a registered trademark of FUJIFILM Corporation.



MEDIVATORS Inc. Headquarters
14605 28th Avenue North
Minneapolis, MN
55447-4822 USA
Tel: +1.800.444.4729
Fax: +1.763.553.3387

MEDIVATORS BV
Sourethweg 11
6422 PC Heerlen
The Netherlands
Tel: +31.45.5.471.471
Fax: +31.45.5.429.695

MEDIVATORS Asia/Pacific Pte. Ltd.
1 A International Business Park
#05-01 Singapore 609933
Tel: +65.6227.9698
Fax: +65.6225.6848

MEDIVATORS Beijing Representative Office
Room 708, 7th Floor Kaiheng Center, Block B
No. 2 Chaoyangmennei Street
Dongcheng District, Beijing China 100010
Tel: +8610.6567.8446
Fax: +8610.6567.8445

www.medivators.com

© MEDIVATORS Inc. P/N 50098-506