


# | Disposable Bougie Dilator

✓ Dilation



## Usage

- It is used for treating with esophageal stenosis and cardia stenosis.

## Features

- Made by medical macromolecule, with proper rigidity.
- Equipped with 2 radiopaque marks on each bougie, easy to observe under X-ray.
- Various bougie sizes, meets different requirements.

## Specifications

Model	Bougie diameter (mm)	Model	Bougie diameter (mm)
VDK-BO-05	5	VDK-BO-13	13
VDK-BO-06	6	VDK-BO-14	14
VDK-BO-07	7	VDK-BO-15	15
VDK-BO-08	8	VDK-BO-16	16
VDK-BO-09	9	VDK-BO-17	17
VDK-BO-10	10	VDK-BO-18	18
VDK-BO-11	11	VDK-BO-19	19
VDK-BO-12	12	VDK-BO-20	20