


# | Disposable Cytology Brush

- Biopsy



## Usage

- This product is used for cell tissue sampling from respiratory tract and alimentary tract.

## Specifications

Model	Catheter O.D. (mm)	Endoscope Channel (mm)	Working Length (mm)	Brush O.D. (mm)	Brush length (mm)	Characteristics
VDK-BC-18-120-A	1.8	≥2.0	1200	2	10	A:Anti- polluted protective tube
VDK-BC-18-105-A	1.8	≥2.0	1050	2	10	
VDK-BC-18-105-B	1.8	≥2.0	1050	2	10	B:"U" Shaped Brush
VDK-BC-18-160-B	1.8	≥2.0	1600	2	10	
VDK-BC-18-120-B	1.8	≥2.0	1200	2	10	
VDK-BC-23-230-B	2.3	≥2.8	2300	2.85	8	
VDK-BC-23-160-B	2.3	≥2.8	1600	2.85	8	
VDK-BC-23-105-B	2.3	≥2.8	1050	2.85	8	
VDK-BC-18-105-C	1.8	≥2.0	1050	3	10	C: Straight Shaped Brush
VDK-BC-18-120-C	1.8	≥2.0	1200	3	10	
VDK-BC-18-180-C	1.8	≥2.0	1800	3	10	
VDK-BC-18-160-C	1.8	≥2.0	1600	3	10	
VDK-BC-18-115-C	1.8	≥2.0	1150	3	10	
VDK-BC-18-110-C	1.8	≥2.0	1100	3	10	
VDK-BC-18-110-C-2	1.8	≥2.0	1100	2	10	
VDK-BC-23-230-C	2.3	≥2.8	2300	3.5	8	
VDK-BC-23-105-C	2.3	≥2.8	1050	3.5	8	
VDK-BC-23-160-C	2.3	≥2.8	1600	3.5	8	
VDK-BC-23-180-C	2.3	≥2.8	1800	3.5	8	