


| Disposable Polyp Snare

ESD

Polyp Retrieval


EMR




cold snare



Crescent



Oval



Hexagonal

***VEDSNARE,
YOUR EMR FIRST CHOICE***

Usage

- This instrument is used in combination with endoscope for cutting polyps or other redundant tissues with high-frequency current in digestive tract.

Features

- Rotatable type optional, guaranteed one-one rotation.
- Moderated soft wire avoids slipping off.

Specifications

Oval Snare

Model	O.D.(mm)	Length(mm)	Loop Width(mm)	Characteristics
VDK-SD-18-160-20-A1	1.8	1600	20	Standard Rotation
VDK-SD-23-160-10-A1	2.3	1600	10	
VDK-SD-23-230-10-A1	2.3	2300	10	
VDK-SD-23-160-15-A1	2.3	1600	15	
VDK-SD-23-230-15-A1	2.3	2300	15	
VDK-SD-23-160-20-A1	2.3	1600	20	
VDK-SD-23-230-20-A1	2.3	2300	20	
VDK-SD-23-160-25-A1	2.3	1600	25	
VDK-SD-23-230-25-A1	2.3	2300	25	
VDK-SD-23-160-30-A1	2.3	1600	30	
VDK-SD-23-230-30-A1	2.3	2300	30	
VDK-SD-23-160-35-A1	2.3	1600	35	
VDK-SD-23-230-35-A1	2.3	2300	35	
VDK-SD-23-160-50-A1	2.3	1600	50	
VDK-SD-23-230-50-A1	2.3	2300	50	
VDK-SD-23-230-10-A1R	2.3	2300	10	Synchronous Rotation
VDK-SD-23-230-15-A1R	2.3	2300	15	
VDK-SD-23-230-20-A1R	2.3	2300	20	
VDK-SD-23-230-30-A1R	2.3	2300	30	

Crescent Snare

Model	O.D.(mm)	Length(mm)	Loop Width(mm)	Characteristics
VDK-SD-18-160-20-B1	1.8	1600	20	Standard Rotation
VDK-SD-23-160-10-B1	2.3	1600	10	
VDK-SD-23-230-10-B1	2.3	2300	10	
VDK-SD-23-160-15-B1	2.3	1600	15	
VDK-SD-23-230-15-B1	2.3	2300	15	
VDK-SD-23-160-20-B1	2.3	1600	20	
VDK-SD-23-230-20-B1	2.3	2300	20	
VDK-SD-23-160-25-B1	2.3	1600	25	
VDK-SD-23-230-25-B1	2.3	2300	25	
VDK-SD-23-160-30-B1	2.3	1600	30	
VDK-SD-23-230-30-B1	2.3	2300	30	

Hexagonal Snare

Model	O.D.(mm)	Length(mm)	Loop Width(mm)	Characteristics
VDK-SD-18-160-20-C1	1.8	1600	20	Standard Rotation
VDK-SD-23-160-10-C1	2.3	1600	10	
VDK-SD-23-230-10-C1	2.3	2300	10	
VDK-SD-23-160-15-C1	2.3	1600	15	
VDK-SD-23-230-15-C1	2.3	2300	15	
VDK-SD-23-160-20-C1	2.3	1600	20	
VDK-SD-23-230-20-C1	2.3	2300	20	
VDK-SD-23-160-25-C1	2.3	1600	25	
VDK-SD-23-230-25-C1	2.3	2300	25	
VDK-SD-23-160-30-C1	2.3	1600	30	
VDK-SD-23-230-30-C1	2.3	2300	30	
VDK-SD-23-230-20-C1R	2.3	2300	20	
VDK-SD-23-230-30-C1R	2.3	2300	30	